

Beaver[®] Cataract Knives

Trusted Performance, Precision and Innovation

Beaver-Visitec International offers a wide range of conventional knives for cataract surgery. Beaver[®] Slit, Sideport, Crescent and Micro-Sharp[®] Knives come in a wide range of sizes to meet the needs of ophthalmic surgeons.



For ordering information call:
1.866.906.8080/1.781.906.8080
or visit www.beaver-visitec.com

bvi Beaver Visitec
Keeping Your Vision in Sight

Beaver® Cataract Knives

A legacy of excellence. Beaver surgical blades have been synonymous with quality, precision and innovation since 1932.

Beaver® Xstar® Slit Knives



Beaver® Xstar® Slit Knives are available in a wide range of sizes designed for smooth penetration and accurate incision. Packaged in a box of 10.



Description	Single Bevel	Double Bevel
1.8 mm, angled 45 degrees	373718	-
2.2 mm, angled 45 degrees	373722	370822
2.4 mm, angled 45 degrees	373724	370824
2.5 mm, angled 45 degrees	373725	370825
2.6 mm, angled 45 degrees	373726	-
2.65 mm, angled 45 degrees	-	370826
2.75 mm, angled 45 degrees	373727	370827
2.8 mm, angled 45 degrees	373728	370828
3.0 mm, angled 45 degrees	373730	370830
3.2 mm, angled 45 degrees	373732	-
3.5 mm, angled 45 degrees	373735	-

Beaver® Sideport/MVR Knives



Beaver® Sideport/MVR Knives are designed to create precise, consistent and fluid tight stab incisions.




Description	Straight	Angled
0.75 mm (23G), 5 per box	585231	585232
1.0 mm, 10 per box	-	373710
1.15 mm (20G), 10 per box	376630	376620
1.4 mm (19G), 5 per box	585240	-

Beaver® Xstar® Crescent Knives



Beaver® Xstar® Crescent Knives offer the flexibility to create optimal tunnel architecture. Packaged in a box of 10.



Description	Item #
2.5 mm, angled 55 degrees, bevel up	373807
2.5 mm, angled 55 degrees, bevel down	373809
2.5 mm, straight, bevel up	373808

Beaver® Micro-Sharp® Blades & Micro-Unitome™ Knives



Stab blades are color coded for simple and accurate identification. Offered in a wide range of sizes. Knives packaged in a box of 10, Blades packaged in a box of 6.



Description	Knives	Blades
15 Degrees, 1.5 mm Depth (Orange)	377512	377511
15 Degrees, 3.0 mm Depth (Blue)	377514	377513
15 Degrees, 5.0 mm Depth (Green)	377516	377515
22.5 Degrees, 4.0 mm Depth (Light Blue)	377523	377522
30 Degrees, 3.5 mm Depth (Black)	377531	377530
45 Degrees, 3.0 mm Depth (Tan)	377546	377545

Beaver® Optimum™ Knives and Blades



Stab blades are color coded for simple and accurate identification. Offered in a wide range of sizes. Knives packaged in a box of 10, Blades packaged in a box of 6.



Description	Knives	Blades
15 Degree Straight (Green)	374891	374881
15 Degree Bent 30 degrees (Green)	374890	
22.5 Degree Straight (Slate blue)	374892	374882
30 Degree Straight (Black)	374893	374883
45 Degree Straight (Tan)	374894	374884

Beaver® Xstar® Implant Knives



Beaver® Xstar® Implant Knives are designed specifically for widening the incision for lens implantation. Packaged in a box of 10.



Description	Item #
3.2 mm, angled 45 degrees, bevel up	373432
3.5 mm, angled 45 degrees, bevel up	373835
4.1 mm, angled 45 degrees, bevel up	373741
5.2 mm, angled 45 degrees, bevel up	373752
3.5 mm, straight	373435
4.1 mm, straight	373341
5.2 mm, straight	373352

Beaver-Visitec International offers a broad range of high-quality precise products for ophthalmic and specialty microsurgical procedures—from knives, cannulae and fluid control products, electrosurgical devices, and microsurgical instruments to drapes and accessories. Our products are sold in shelf packs and many can be combined in custom procedure packs or instrument trays for added convenience and efficiency.



Beaver-Visitec International, Inc.

411 Waverley Oaks Road, Waltham, MA 02452

Tel: 1.866.906.8080 Fax: 1.866.906.4304

www.beaver-visitec.com

BVI, BVI Logo and all other trademarks (unless noted otherwise) are property of a Beaver-Visitec International ("BVI") company © 2013 BVI