

# Beaver<sup>®</sup> Otolaryngology Products

Trusted Precision for Specialized Procedures

Beaver-Visitec International offers a wide range of Beaver<sup>®</sup> Ear, Nose and Throat (ENT) blades designed to help enhance control and operative efficiency in a wide range of procedures.

ENT/Head & Neck Procedure Blades  
Myringotomy Blades and Knives  
Tympanoplasty/Stapedectomy Blades  
Custom ENT Procedure Kits/Packs



For ordering information call:  
1.866.906.8080/1.781.906.8080  
or visit [www.beaver-visitec.com](http://www.beaver-visitec.com)








**bvi** Beaver Visitec  
*Keeping Your Vision in Sight*

# Beaver® Otolaryngology Products


## ENT/Head & Neck Procedure Blades

Beaver® ENT/Head & Neck blades are designed to help enhance control and operative efficiency in a wide range of procedures, from facial nerve dissection to tonsillectomy. Because the requirements of each product are so unique and specialized, each Beaver® blade is shaped, angled and honed to help achieve the desired clinical result.

### Beaver® Mini-Blade®

Description	Item #
 sharp on two sides, pointed tip, 6 per box	376300
 rounded tip, sharp on one side, 24 per box	376400
 rounded narrow tip, sharp on one side, 24 per box	376401
 smaller alternative to #11 general blade, 24 per box	376500
 smaller alternative to #15 general blade, 24 per box	376700
 straight, sharp all around, 180° cutting surface double bevel blade, 6 per box	376900
 angled 60°, sharp all around, 180° cutting surface double bevel blade, 6 per box	376910

### Beaver® Mini-Blade® Knife

Description	Item #
 straight, sharp all around, 180° cutting surface double bevel blade, 3 per box	374769



### Tonsillectomy Blade

Description	Item #
 sickle sharp, 24 per box	376800

### Needle Blade

Description	Item #
 1.5 mm cutting edge, 3 per box	375940
 3.0 mm cutting edge, 3 per box	375910
 5.0 mm cutting edge, 3 per box	375920
 7.0 mm cutting edge, 3 per box	375930

### Beaver® Handles

Description	Item #
 Round, Knurled, Stainless 10.0 cm L x .50 cm D, 1 per box	374310
 Round, Knurled, Stainless 13.0 cm L x .50 cm D, 1 per box	371350

\* Beaver Handles will accommodate all Beaver blades.




## Myringotomy Blades and Knives

Beaver® Myringotomy blades and knives are intended to offer superior precision and control in delicate, typically pediatric myringotomy procedures. A narrower shaft has been engineered to improve visibility of the surgical field for greater intraoperative control. Tips designed for right or left cutting, or two-directional cutting, provide the versatility needed for fine tympanic incisions.

Myringotomy Blades	
Description	Item #
 spear tip, offset 45°, narrow shaft, 6 per box	377120
 lance tip, downward cut, offset 45°, 6 per box	377100
 lance tip, upward cut, offset 45°, 6 per box	377110
 lance tip, downward cut, narrow shaft, offset 45°, 6 per box	377121
Myringotomy Knives	
Description	Item #
 spear tip, offset 45°, narrow shaft, 3 per box	377122
 lance tip, downward cut, narrow shaft, offset 45°, 3 per box	377124

## Tympanoplasty/Stapedectomy Blades

A range of Beaver® blade designs, shapes and angles is offered for the special requirements of both tympanoplasty and stapedectomy, from general canal procedures to elevating the delicate tympanic membrane. A narrower shaft has been engineered to enhance visibility of the operative site.

Tympanoplasty Blades	
Description	Item #
 60° angled, sharp all around 2.5 mm W, bevel down, 3 per box	377200
 50° angled, sharp all around 2.0 mm W, bevel down, 3 per box	374451
 straight, sharp all around 2.5 mm W, , 3 per box	377210
 Sickle Edge Blade, straight 84 mm L 5.0 mm cutting surface, 3 per box	377300

